# **BIDS AND AWARDS COMMITTEE**

Website: www.rsu.edu.ph

Tel No. (042) 567 - 5273/5859/6234

Tele Fax No. (042) 567-5270

REQUEST FOR QUOTATION Solicitation No. RSU-21-18 Quotation No. 18-06-054

Date	: JUNE 26, 2018
Company Name: Address:	
Audress	

Sir/Madame:

Please quote your lowest price on the Item/s below, subject to the General Conditions and Eligibility Requirements for the Suppliers, stating the shortest time of the delivery and submit your quotation duly signed by your representative not later than <u>JUNE 29, 2018</u> at exactly 5:00 in the afternoon in the return envelope attached herewith.

**Project**: Procurement of Various Laboratory Supplies of Science Laboratory

**ABC:** PhP143,267.00

**QTY:** 8730

L Acetone 2	UNIT	ITEM DESCRIPTION	QTY	ESTIMATED UNIT COST	ESTIMATED COST
kg         Alum. Powder         1           L         Barfoed's solution         1           vial         BS A standard (bovine serum albumin)         1           L         Buffer ph-standard 4         1           L         Buffer ph-standard 7         1           g         Calcium chloride         500           vial         Capillary tube         2           pc         Cly triangle         20           pc         Cord, rassorted size (#9,10,16)         100           pc         Cover slip         50           pc         Cover slip         50           pc         Cover slip         50           pc         Covible and cover, 30mL         12           L         Dimethyl sulfoxide         2.5           L         Ethanol, 58 technical grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           L         Hexane         1           L         Hydrochloric acid         2.5           L         Hydrochloric acid         2.5           L         Hydrochloric acid	L	Acetone	2		
L   Barfoed's solution	g	Agar-Agar	500		
vial         BSA standard (bovine serum albumin)         1           L         Buffer ph-standard 10         1           L         Buffer ph-standard 4         1           L         Buffer ph-standard 7         1           g         Calcium chloride         500           vial         Capillary tube         2           pc         Clay triangle         20           pc         Core, assorted size (#9,10,16)         100           pc         Cover slip         50           pc         Crucible and cover, 30mL         12           L         Dimethyl sulfoxide         2.5           L         Ethanol, AR grade         2.5           L         Ethanol, Sey technical grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           L         Hydrochloric acid         2.5           L         Hydrochloric acid         2.5           L         Hydrochloric acid         1           L         Hydrochloric acid         1           L         Hydrochloric acid         2.5           L         Hy	kg	Alum. Powder	1		
L Buffer ph-standard 10	L	Barfoed's solution	1		
L Buffer ph-standard 4	vial	BSA standard (bovine serum albumin)	1		
L Buffer ph-standard 7	L	Buffer ph-standard 10	1		
Section   Calcium chloride   Soo   Soo	L	Buffer ph-standard 4	1		
Section   Calcium chloride   Soo   Soo	L	Buffer ph-standard 7	1		
vial         Capillary tube         2           pc         Clay triangle         20           pc         Cork, assorted size (#9,10,16)         100           pc         Cover slip         50           pc         Crucible and cover, 30mL         12           L         Dimethyl sulfoxide         2.5           L         Ethanol, AR grade         2.5           L         Ethanol, AR grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Marconkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g <td>g</td> <td></td> <td>500</td> <td></td> <td></td>	g		500		
pc         Clay triangle         20           pc         Cork, assorted size (#9,10,16)         100           pc         Cover slip         50           pc         Crucible and cover, 30mL         12           L         Dimethyl sulfoxide         2.5           L         Ethanol, AR grade         2.5           L         Ethanol, 95% technical grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrochloric acid         2.5           L         Hydrochloric acid         2.5           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           g         Mac conkey agar         500           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           p <td< td=""><td></td><td>Capillary tube</td><td>2</td><td></td><td></td></td<>		Capillary tube	2		
pc         Cork, assorted size (#9,10,16)         100           pc         Cover slip         50           pc         Crucible and cover, 30mL         12           L         Dimethyl sulfoxide         2.5           L         Ethanol, AR grade         2.5           L         Ethanol, 95% technical grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g Iron(II) sulfate         500           g Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Manitol Salt agar         500           g         Mariand standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Minhydrin solution         1           L         Minhydrin sol	рс		20		
pc         Cover slip         50           pc         Crucible and cover, 30mL         12           L         Dimethyl sulfoxide         2.5           L         Ethanol, AR grade         2.5           L         Ethanol, 95% technical grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrogen peroxide         1           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Mac conkey agar         500           g         Manitol Salt agar         500           g         Maridol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g	рс		100		
pc         Crucible and cover, 30mL         12           L         Dimethyl sulfoxide         2.5           L         Ethanol, AR grade         2.5           L         Ethanol, 95% technical grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           d         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500	рс		50		
L         Dimethyl sulfoxide         2.5           L         Ethanol, AR grade         2.5           L         Ethanol, 95% technical grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Mac conkey agar         500           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L	рс	•			
L       Ethanol, AR grade       2.5         L       Ethanol,95% technical grade       6         pc       Glass tubing, 7mm x 5 ft.       10         ml       Gram staining       250         L       Hexane       1         tube       High vacuum grease (150g), Silver grease       1         L       Hydrochloric acid       2.5         L       Hydrogen peroxide       1         g       Iron(II) sulfate       500         gnack       Lens cleaning tissues, 100x150mm (100sheets)       2         g       Mac conkey agar       500         g       Manitol Salt agar       500         set       Mcfarland standards       1         pc       Medicine dropper       20         g       Metal zinc       500         g       Metal zinc       500         L       Methyl alcohol AR       2.5         g       Mueller Hinton Agar       500         L       Ninhydrin solution       1         g       Nutrient agar       500         g       Nutrient broth       500         g       Petri dish (glass) 15 x 100mm,       50         g       Potassium chlorate <td></td> <td></td> <td></td> <td></td> <td></td>					
L         Ethanol,95% technical grade         6           pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           g         Iron(II) sulfate         500           g         Mac conkey agar         500           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient broth         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm, </td <td></td> <td></td> <td></td> <td></td> <td></td>					
pc         Glass tubing, 7mm x 5 ft.         10           ml         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           g         Iron(II) sulfate         500           g         Mac conkey agar         500           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         2.5           L         Ninhydrin solution         1           g         Nutrient broth         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium Chlorate	L				
mI         Gram staining         250           L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500	рс				
L         Hexane         1           tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Macin Jala agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500					
tube         High vacuum grease (150g), Silver grease         1           L         Hydrochloric acid         2.5           L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500           kg         Potassium Chromate         500					
L       Hydrochloric acid       2.5         L       Hydrogen peroxide       1         g       Iron(II) sulfate       500         pack       Lens cleaning tissues, 100x150mm (100sheets)       2         g       Mac conkey agar       500         g       Manitol Salt agar       500         set       Mcfarland standards       1         pc       Medicine dropper       20         g       Metal zinc       500         L       Methyl alcohol AR       2.5         g       Mueller Hinton Agar       500         L       Ninhydrin solution       1         g       Nutrient agar       500         g       Nutrient broth       500         roll       Parafilm       1         pc       Petri dish (glass) 15 x 100mm,       50         g       Potassium chlorate       500         kg       Potassium Chromate       500	tube				
L         Hydrogen peroxide         1           g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500           kg         Potassium Chromate         500					
g         Iron(II) sulfate         500           pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500           kg         Potassium Chromate         500	L				
pack         Lens cleaning tissues, 100x150mm (100sheets)         2           g         Mac conkey agar         500           g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500           kg         Potassium Chromate         500					
g       Mac conkey agar       500         g       Manitol Salt agar       500         set       Mcfarland standards       1         pc       Medicine dropper       20         g       Metal zinc       500         L       Methyl alcohol AR       2.5         g       Mueller Hinton Agar       500         L       Ninhydrin solution       1         g       Nutrient agar       500         g       Nutrient broth       500         roll       Parafilm       1         pc       Petri dish (glass) 15 x 100mm,       50         g       Potassium chlorate       500         kg       Potassium Chromate       500					
g         Manitol Salt agar         500           set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500           kg         Potassium Chromate         500	•				
set         Mcfarland standards         1           pc         Medicine dropper         20           g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500           kg         Potassium Chromate         500					
pc         Medicine dropper         20					
g         Metal zinc         500           L         Methyl alcohol AR         2.5           g         Mueller Hinton Agar         500           L         Ninhydrin solution         1           g         Nutrient agar         500           g         Nutrient broth         500           roll         Parafilm         1           pc         Petri dish (glass) 15 x 100mm,         50           g         Potassium chlorate         500           kg         Potassium Chromate         500					
L Methyl alcohol AR  g Mueller Hinton Agar  L Ninhydrin solution  g Nutrient agar  g Nutrient broth  roll Parafilm  pc Petri dish (glass) 15 x 100mm,  g Potassium chlorate  kg Potassium Chromate  2.5  500  1  1  500  500  600  600  600  6	· ·				
gMueller Hinton Agar500LNinhydrin solution1gNutrient agar500gNutrient broth500rollParafilm1pcPetri dish (glass) 15 x 100mm,50gPotassium chlorate500kgPotassium Chromate500					
L       Ninhydrin solution       1       1         g       Nutrient agar       500       500         g       Nutrient broth       500       500         roll       Parafilm       1       1         pc       Petri dish (glass) 15 x 100mm,       50       50         g       Potassium chlorate       500       500         kg       Potassium Chromate       500       500					
g         Nutrient agar         500	L				
gNutrient broth500rollParafilm1pcPetri dish (glass) 15 x 100mm,50gPotassium chlorate500kgPotassium Chromate500	g				
rollParafilm1pcPetri dish (glass) 15 x 100mm,50gPotassium chlorate500kgPotassium Chromate500					
pcPetri dish (glass) 15 x 100mm,50gPotassium chlorate500kgPotassium Chromate500					
g Potassium chlorate 500 kg Potassium Chromate 500					
kg Potassium Chromate 500					
. R T FOLGSSIUM JOULE I DUO I I	g	Potassium iodide	500		

SVP (Section 53.9)



## **BIDS AND AWARDS COMMITTEE**

Website: www.rsu.edu.ph

Tel No. (042) 567 - 5273/5859/6234

Tele Fax No. (042) 567-5270

g	Potato dextrose agar		500		1
roll	Rubber Tubing		2		
g	Sodium Chloride		500		
g	Sodium citrate, Trisodium citrate		500		
g	Sodium Hydroxide		500		
L	Sulfuric acid		1		
рс	Thermometer, alcohol based, 200 oC		20		
рс	Watch glass 100mm		20		
рс	Watch glass 75mm		20		
рс	Wire gauze ceramic center		120		
L	Wright stain		1		_
TOTAL		QTY	8730	ESTIMATED COST	

Υ	O	u	rs	

Prof. MARIO A. FETALVER, Jr., Ph. D.

Chairperson, BAC	
Canvassed by:	
<b>LARRY I. FIRMALO</b> Procurement Officer	
BIDS & AWARDS COM ROMBLON STATE UNI Odiongan, Romblon	• •
Sir/Madame:	
Delivery Period Warranty Price Validity	: : :
After having carefully i	read accepted your General Conditions, I/We quote you on the item at prices noted above.
Printed Name / Signat	<del></del> ure
Tel. No. /Cellphone No Email Address	o
Dutc	

### **REQUEST FOR QUOTATION**

### **Terms & Conditions:**

- 1. All entries must be type/hand-written.
- 2. Delivery period within 15 calendar days.
- 3. Warranty shall be for a period of six (6) months for supplies and materials, one (1) year for equipment, from the date of acceptance by the procuring entity.
- 4. Price validity shall be for a period of three months.





# **BIDS AND AWARDS COMMITTEE**

Website: www.rsu.edu.ph

Tel No. (042) 567 - 5273/5859/6234

Tele Fax No. (042) 567-5270

- 5. PhilGEPS registration Certificate shall be attached upon submission of the quotation.
- 6. Bidder shall submit original brochures showing certifications of the production being offered.
- 7. Payment shall be made after the inspection.
- 8. Deliver items at Supply Office, Romblon State University, Odiongan, Romblon (042-567-5375)

#### **CHECKLIST OF ELIGIBILITY REQUIREMENTS FOR THE BIDDERS/SUPPLIERS:**

The Eligibility Envelope shall contain the following:

- 1. Mayor's/Business permit
- 2. PhilGEPS Registration Number
- 3. Income/Business Tax Return
- 4. Omnibus Sworn Statement

Note: Submitted documents must be properly authenticated.